二重 β 崩壊検出のためのCANDLES の開発

Kayoko Ichihara

Department of Physics, Osaka Univ.

T.Kishimoto, I.Ogawa, R.Hazama, S.Yoshida, S.Umehara, K.Matsuoka and Y.Hirano

for the CANDLES collaborations

Contents

Physics motivation of double beta decay
 Study of double beta decay in the world

3. Double beta decay of ⁴⁸Ca

- 4. CANDLES system
- 5. BG rejection (Pulse shape discrimination)

6. Goal of CANDLES system

Physics motivation of double beta decay - 1 -

1-1. Introduction

Neutrino oscillation \longrightarrow mixing angle, Δm^2

- Absolute mass of neutrino? & Mass pattern?
- Why neutrinos are so light ?
- Is there heavy neutrino ? Dirac or Majorana ?
- Matter-antimatter asymmetry in Universe.

Is lepton number violated ?

Most sensitive prove

New physics beyond the Standard Model !!

1-2. Double beta decay



1-3. Double beta decay

⁴⁸Ca,⁷⁶Ge,⁸⁷Se,⁹⁶Zr,¹⁰⁰Mo, ¹¹⁶Cd, ¹²⁸Te,¹³⁰Te,¹³⁶Xe,¹⁵⁰Nd.





1-4. Effective mass vs minimum mass



http://www.aps.org/neutrino

Double beta decay in the world - 2 -

2-1. Current situation in experimental DBD

Isotope	Exposure	BG	$T0^{\nu\beta\beta}_{1/2}$	<m<sub>\beta\beta\beta}></m<sub>	Collaboration
	(kmole-y)	(counts)	(year)	(eV)	(Inclusive/ Exclusive)
⁴⁸ Ca	5x10 ⁻⁵	0	> 1.4 x 10 ²²	< 7.2-44.7	ELEGANT VI
⁷⁶ Ge	0.943	61	$= 1.2 \text{ x } 10^{25}$	= 0.44	Helberg-Moscow
⁷⁶ Ge	0.117	3.5	$> 1.6 \text{ x } 10^{25}$	< 0.33-1.35	1 EEX
⁸² Se	7x10 ⁻⁵	0	$> 2.7 \text{ x } 10^{22} \text{ *}$	< 5.0	(HC) *: 68%CL
¹⁰⁰ Mo	5x10-4	4	$> 5.5 \text{ x } 10^{22}$	< 2.1	E GANT V
¹¹⁶ Cd	1x10 ⁻³	14	$> 1.7 \text{ x } 10^{23}$	< 1.7	(Ex)
¹²⁸ Te	Geochem.	NA	$> 7.7 \text{ x } 10^{24}$	< 1.1-1.5	(In)
¹³⁰ Te	0.0025	5	$> 5.5 \text{ x } 10^{23}$	< 0.37-1.9	CUORICINO (In)
¹³⁶ Xe	7x10 ⁻³	16	$> 4.4 \text{ x } 10^{23}$	< 1.8-5.2	Gotthard
¹⁵⁰ Nd	6x10 ⁻⁵	0	$> 1.2 \text{ x } 10^{21}$	< 3.0	(Fer

Uncertainty in $|M_{0v}|^{(Ex)}$

http://www.aps.org/neutrino

Z-Z. Heiderberg-woscow (Enriched "Ge



Phys. Lett. B 586 (2004) 198

2-3. CUORICINO (Cryogenic Bolometor)



 $T_{1/2}^{0\nu} > 5.5x10^{23}years$ for ¹³⁰Te(~30kg) <mv> < 0.37 - 1.9eV Phys. Lett. B 584 (2004) 260



2-5. Present status & future experiment



Double beta decay of ⁴⁸Ca - 3 -

3-1. Double beta decay of ⁴⁸Ca

- CaF₂ scintillation detector Source= Detector
 - : High Detection Efficiency
- Largest Q_{ββ} value : Low background



3-2. ELEGANT VI



- Passive shield for neutron
 - : LiH+paraffin 15mm
 - : Cd sheet 0.6mm
 - : $H_3BO_3+H_2O$ tank

<u>ELE</u>ctron <u>GA</u>mma ray <u>N</u>eutrino <u>T</u>elescope

Detector & Source
 : CaF₂ (Eu)
 45x45x45 mm³, 23crystals 6.7kg

4π active shield

: CaF₂(pure):active L.G.

45x45x200 mm³, 48crystals

- **: Csl(Tl)** 65x65x250mm³, 38modules
- Passive shield for γ ray
 - : OFHC Cu 5cm
 - : Pb 10cm

3-3. Oto Cosmo Observatory



3-4. BG free measurement @ ELE-VI 100 counts/20 keV **Energy resolution** : 3.1(FWHM) @4.27MeV 10² **Detection eff.** : 58% from simulation 10 $T_{1/2}^{0v} > 1.4 \times 10^{22} y$ 1 ⁴⁸Ca: 4271 16C 36 30 104000 5000 2000 3000 Energy (keV) NPA730(2004)215

CANDLES system - 4 -

4-1. Improvement of sensitivity

ELEGANT VI (6.6g of ⁴⁸Ca)@Oto
 4π active shield + Largest Q_{ββ} value

BG free measurement

Increase the number of ⁴⁸Ca nuclei (Order of kg)
 Large volume detector

 CANDLES

<u>CA</u>lcium fluoride for studies of <u>N</u>eutorino

and <u>D</u>ark matter by <u>L</u>ow <u>E</u>nergy <u>S</u>pectrometer

• 4π active shield + largest $Q_{\beta\beta}$ value

Huge volume detector !!

4-2. CANDLES system

<u>CA</u>lcium fluoride for studies of <u>N</u>eutorino and <u>D</u>ark matter by <u>L</u>ow <u>E</u>nergy <u>S</u>pectrometer

Undoped CaF₂ scintillator

: Long attenuation length

- Liquid scintillator
 - : 4π active shield

Decay time of signal 900nsec : CaF₂(pure)

A few 10nsec : Liq. Scinti.





4-4. Progress of CANDLES

CANDLES I

Light collection efficiency : Optimization of mixture

BG rejection : Pulse shape discrimination Liq./CaF₂, α/γ , double pulse

CANDLES III 60crystals x 3.2kg = 191kg

Under Construction

R&D study

- Reduction of radioactive impurities in CaF₂ crystal
- Enrichment ⁴⁸Ca

BG rejection Pulse shape discrimination - 5 -

5-1. Background







 α/γ & double pulse rejection

5-2. Set Up (α/γ discrimination)



5-3. Reference Pulse





Goal of CANDLES system - 6 -

6-1. Sensitivities of CANDLES

CANDLES III ————



< mv > : 0.11 eV

- CaF_2 : 60 crystals x 3.2 kg 191 kg
- E-Resolution: 4%@4.27MeV
- Measurement time: 5year
- CANDLES IV
 - CaF₂: 1000crystals x 3.2kg 3.2ton
 - E-Resolution: 3.5%@4.27MeV
 - Measurement time: 10year
- CANDLES V
 - Enrichment of ⁴⁸Ca 0.187% => 2%

or Scale up detector (40ton) & BG reduction

Measurement time: 10year

Summary

0νββ decay

Lepton number violation Majorana absolute mass of neutrino



- BG free measurement @ ELE-VI successful
 Increase the number of ⁴⁸Ca nuclei CANDLES
- Performance of CANDLES system
 - Light collection efficiency : good
 - Pulse shape discrimination Liq./CaF₂, double pulse, α/γ : good
- CANDLES III (Sensitivity <m_n> : 0.56eV) Under Construction
 Future prospect CANDLES V ----> <m_n> : 0.04eV

Thank you for your attention !!